



GotchaDwi.com

Lt. Allan Marx

479-629-4986

amarx@premierbiotech.com

P.O. Box 180074

Fort Smith, AR 72918

You Can't Beat Our Oral Fluid Drug Test



No products are available that will allow a donor to cheat this test

ADVANTAGES OF ORAL FLUID TESTING:

- Eliminates the need for controlled collection sites
- Prevents adulteration by substitution and dilution

ORAL FLUID TESTING

Only requirement: Donor must have nothing in the mouth 10 minutes prior to collection.

URINE

Substitution

- Donor-obtained urine
- Synthetic urine or other substitute substance

Internal Dilution

- Drinking large amounts of liquids
- Golden Seal Tea

External Diltuion

- Adding water to urine sample

** If enough water is added internally or externally, the same may test negative.*

Products like these below do not exist for oral fluid testing:



SIX PANEL

AMP
amphetamine
50 ng/ml

COC
cocaine
20 ng/ml

MET
methamphetamine
50 ng/ml

OPI
opiates
40 ng/ml

PCP
phencyclidine
10 ng/ml

THC
marijuana
50 ng/ml

TEN PANEL DEVICE ADDS

BAR
barbiturates
50 ng/ml

BZO
benzodiazepines
10 ng/ml

MTD
methadone
30 ng/ml

PPX
propoxyphene
50 ng/ml